

BASIS FLOW II

- Multi-Purpose Pourable Cold Curing Acrylic Denture Resin
- Cures in Pressure Pot in 30-60 Minutes at 131°F
- Long Working Time
- Japanese Technology



BASIS FLOW II



www.ydp-usa.com
www.yamahachidentalusa.com
Phone: (877) 320-2909
info@yamahachidentalusa.com

yamahachi 
Dental Products USA

BASIS FLOW II is a Brand distributed by Yamahachi Dental Products USA Inc.







BASIS FLOW II MAIN FEATURES

- Multi-Purpose Pourable Cold Curing Acrylic Denture Resin
- Cures in Pressure Pot in 30-60 Minutes at 131°F
- Long Working Time
- Japanese Technology



BASIS FLOW II TECHNICAL FEATURES

CURING	
Powder / Liquid mixing ratio	100 mg / 60 ml
Pressure Pot	55°C / 29 PSI
PHYSICAL PROPERTIES	
Flexural Strength	90 Mpa
Elasticity	14.6 mm
Vickers Hardness	15.6 Hv
SHADES	
O Pink	
V Pink	
LF Pink	
Clear	

Buy BASIS FLOW II

<https://yamahachidentalusa.com/products/basis-flow-ii-acrylic>

Package information:

Powder: 500 g

Liquid: 500 ml

BASIS FLOW II



www.ydp-usa.com
www.yamahachidentalusa.com
Phone: (877) 320-2909
info@yamahachidentalusa.com

yamahachi
Dental Products USA

BASIS FLOW II is a Brand distributed by Yamahachi Dental Products USA Inc.